




Specifications

- Sample Size: 100 µL.
- Sample Type : Whole Blood, Serum, Plasma, Urine etc.
- Sample Application – Syringe, Sample Cup, Collection Tubes, Capillary.
- Analysis Time : 55 Seconds.
- Temperature : Room Temperature, 15°C–32°C (60°F–90°F).
- Dimensions H X W X D : 38.5 X 32.5 X 30.5 Cm.
- Weight : 6 Kgs (13 Lbs) Approx.
- Printout : Built-in, Thermal Roll Printer; 32 Characters Wide
- Electronic : Microprocessor Controlled.
- Relative Humidity : < 85%, Non-condensing.
- Normal Value : Flagging Of Abnormal Results; User Programmable Ranges.
- Calibration : Fully Automatic, 1-point With Every Sample; 2-point Every 4 Hours.
- Power Requirements : 110–240 V, 50/60 Hz (self-adjusting) 1.4 Amp Max.

Features

- Reliable and User Friendly menu.
- Precise and Accurate results.
- Long life, maintenance free electrodes.
- Inbuilt Cleaner (Optional).
- Smart Reagent Kit reactive with electronic chip.
- Extremely Low cost per test.
- Inbuilt Parameter conversion functionality.

Available Parameter

- Sodium | Potassium | Chloride | Calcium (Ion) | Lithium | pH
- Sodium | Potassium | Chloride | Calcium (Ion) | Lithium
- Sodium | Potassium | Chloride | Calcium (Ion) | pH
- Sodium | Potassium | Chloride | Calcium (Ion)
- Sodium | Potassium | Chloride | Lithium
- Sodium | Potassium | Calcium (Ion) | Lithium
- Sodium | Potassium | Chloride | pH
- Sodium | Potassium | Calcium (Ion) | pH
- Sodium | Potassium | Chloride
- Sodium | Potassium | Calcium (Ion)
- Sodium | Potassium



Marketing By



Off: 9987750054

Sai Prasad Society, Shop No. 3, Road No. 2,
Opposite Basera Studio, Singh Estate, Kandivali East,
Mumbai 400101

Manufactured by:

Promea Therapeutics Pvt Ltd.

Plot No. M4, Medical Device Park, Sultanpur (V), Patancheru (M),
Sangareddy (Dist), Telangana - 502319. www.promea.in